




Recombinant Human T cell receptor beta constant 1 (TRBC1), partial

Product Code	CSB-EP024291HU1
Abbreviation	Recombinant Human TRBC1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01850
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	DLNKVFPPEVAVFEPSEAEISHTQKATLVCLATGFFPDHVELSWVWVNGKEVHS GVSTDPQPLKEQPALNDSRYCLSSRLRVSATFWQNPRNHFRQCQVQFYGLSE NDEWTQDRAKPVTQIVSAEAWGRAD
Research Area	Others
Source	E.coli
Target Names	TRBC1
Expression Region	1-129aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	19.7 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.