




# Recombinant Human Endogenous retrovirus group K member 113 Gag polyprotein (HERVK\_113)

<b>Product Code</b>	CSB-EP352942HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P62684
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	GQTKSKIKSKYASYLSFIKILLKRGGVK VSTKNLIKLFQII EQFCPWFPEQGTLDL KDWKRIGKELKQAGRKGNIIPLTVWNDWAIKAALEPFQTEEDSVSVSDAPGSC IIDCNEKTRKKSQKETESLHCEYVAEPVMAQSTQNADYNQLQEVIYPETLKLEG KGPELMGPSESKPRGTSPLPAGQVPVTLQPQKQVKENKTQPPVAYQYWPPA ELQYQPPPE SQYGYPGMPPAPQGRAPYPQPPTRRNLNPTAPPSRQGSSELHEII DKSRKEGDTEAWQFPVTLELMPPGEGAQEGEPPTVEARYKSFSIKMLKDMKE GVKQYGPNSPYMRTL LDSIAHGHRLIPYDWEILAKSSLSPSQFLQFKTWWIDG VQEQVRRNRAANPPVNIDADQLLGIGQNWSTISQQALMQNEAIEQVRAICLRA WEKIQDPGSTCPSFNTVRQGSKEPYPDFVARLQDVAQKSIADKARKVIVELM AYENANPECQSAIKPLKGKVPAGSDVISEYVKACDGMGGAMHKAMLMAQAIT GVVLGGQVRTFGGKCYNCGQIGHLKKNC PVLNKQINITIQTATTTGREPPDLCPR CKKGKHWASQCRSKFDKNGQPLSGNEQRGQPQAPQQTGAFPIQPFVPPQGF QGQQPPLSQVFQGISQLPQYNNCPPPQAAVQQ
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	HERVK_113
<b>Protein Names</b>	Recommended name: HERV-K_19p13.11 provirus ancestral Gag polyprotein Short name= Gag polyprotein Alternative name(s): HERV-K113 Gag protein Cleaved into the following 3 chains: 1. Matrix protein 2. Capsid protein 3. Nucleocapsid protein
<b>Expression Region</b>	2-666aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	81.3 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.