




Recombinant Daucus carota Phytosulfokine-alpha

Product Code	CSB-EP348219DIR
Abbreviation	Recombinant Daucus carota Phytosulfokine-alpha protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P58261
Product Type	Recombinant Proteins
Immunogen Species	Daucus carota(Wild carrot)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	YIYTQ
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Phytosulfokine-alpha Short name= PSK-alpha Cleaved into the following chain: 1. Phytosulfokine-beta Short name= 2. PSK-beta
Expression Region	1-5aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	16.0 kDa
Protein Length	Full length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.