




Recombinant Mouse Histone H2B type 1-A (Hist1h2ba)

Product Code	CSB-EP010401MO
Abbreviation	Recombinant Mouse Hist1h2ba protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P70696
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PEVAVKGATISKKGFKKAVTKTQKKEGRKRKRRCRKESYSIYIKVLKQVHPDTG ISSKAMSIMNSFVTDIFERIASEASRLAHYNKRSTITSREIQTAVRLLLPGELAKH AVSEGTKAVTKYTSSK
Research Area	Others
Source	E.coli
Target Names	Hist1h2ba
Protein Names	Recommended name: Histone H2B type 1-A Alternative name(s): Histone H2B, testis Testis-specific histone H2B
Expression Region	2-127aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	21.5 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.