




Recombinant Lake Victoria marburgvirus Nucleoprotein (NP), partial

Product Code	CSB-EP333984LBR
Abbreviation	Recombinant Lake Victoria marburgvirus NP protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35263
Product Type	Recombinant Proteins
Immunogen Species	Lake Victoria marburgvirus(strain Popp-67)(MARV)(Marburg virus (strain West Germany/Popp/1967))
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	NVGEQYQQLREAAHDAEIKLQRRHEHQEIQAIAEDDEERKILEQFHLQKTEITH SQTAVLSQKREKLARLAAEIENNIVEDQGFKQSQNQVSQSFLNDPTPVEVTV QARPINRPTALPPPVDNKIEHESTEDSSSSSSFVDLNDPFPALLNEDEDTLDDSV MIPSTTSREFQGIPAPPRQSQDLNNSQGKQDESTNPIKKQFLRYQELPPVQE DDESEYTTDSQESIDQPGSDNEQGVDLPPPLYAQEKRPDIQHPAVSSQDP FGSIGDVNGDILEPIRSPSSPSAPQEDTRAREAYELSPDFTNYEDNQQNWPQR VVTKKGRFTFLYPNDLLQTNPPESLITALVEEYQNPVSAKELQADWPDMSFDER RHVAMNL
Research Area	Signal Transduction
Source	E.coli
Target Names	NP
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	317-695aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	50.4 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.