




Recombinant Human Cis-aconitate decarboxylase (ACOD1)

Product Code	CSB-EP011825HU
Abbreviation	Recombinant Human ACOD1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NK06
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MMLKSITESFATAIHGLKVGHLTDRVIQRSKRMILDTLGAGFLGTTTEVFHIASQ YSKIYSSNISSTVWGQPDIRLPPTYAAAFVNGVAIHSMDFDDTWHPATHPSGAV LPVLTALAEALPRSPKFSGLDLLAFNVGIEVQGRLLHFAKEANDMPKRFHPPS VVGTLGSAAAASKFLGLSSTKCREALAIAVSHAGAPMANAATQTKPLHIGNAAK HGIEAAFLAMLGLQGKNQVLDLEAGFGAFYANYSPKVLPSIASYSWLLDQQDV AFKRFPAPHLSTHWWADAAASVRKHLVAERALLPTDYIKRIVLRIPNVQYVNRPF PVSEHEARHSFQYVACAMLLDGGITVPSFHECQINRPQVRELLSKVELEYPPD NLPSFNILYCEISVTLKDGATFTDRSDTFYGHWRKPLSQEDLEEKFRANASKML SWDTVESLIKIVKNLEDLEDCSVLTLLKGPSPEVASNSPACNNSITNLS
Research Area	Others
Source	E.coli
Target Names	ACOD1
Protein Names	Recommended name: Immune-responsive gene 1 protein homolog
Expression Region	1-481aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	58.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.