



Recombinant Human Paired box protein Pax-5 (PAX5)

Product Code	CSB-EP2065HU
Abbreviation	Recombinant Human PAX5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q02548
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MDLEKNYPTPRTSRTGHGGVNQLGGVVFVNGRPLPDVVRQRIVELAHQGVRP CDISRQLRVSHGCVSKILGRYYETGSIKPGVIGGSKPKVATPKVVEKIAEYKRQ NPTMFAWEIRDRLLAERVCDNDTVPSVSSINRIIRTKVQQPPNQPVPASSHSIV STGSVTQVSSVSTDSAGSSYSISGILGITSPSADTNKRKRDEGIQESVPVNGHS LPGRDFLRKQMRGDLFTQQQLEVLDRVFERQHYSIDIFTTTEPIKPEQTTEYSA MASLAGGLDDMKANLASPTPADIGSSVPGPQSYPIVTGRDLASTTLPGYPPHV PPAGQGSYSAPTLTGMVPGSEFSGSPYSHQPQYSSYNDSWRFPNPGLLGSPY YYSAAARGAAPPAAATAYDRH
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	PAX5
Expression Region	1-391aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.