




Recombinant Olea europaea Pollen allergen Ole e 7, partial

Product Code	CSB-EP305742OEC
Abbreviation	Recombinant Olea europaea Pollen allergen Ole e 7 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P81430
Product Type	Recombinant Proteins
Immunogen Species	Olea europaea(Common olive)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	APSQSTVTALLTSCVSYIDDQ
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Pollen allergen Ole e 7Alternative name(s): Allergen Ole e VII Allergen= Ole e 7
Expression Region	1-21aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	17.5 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.