




Recombinant *Klebsiella pneumoniae* subsp. *pneumoniae* UDP-3-O-acyl-N-acetylglucosamine deacetylase (IpxC)

Product Code	CSB-EP408442KAX
Abbreviation	Recombinant <i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> IpxC protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6T4N9
Product Type	Recombinant Proteins
Immunogen Species	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> (strain ATCC 700721 / MGH 78578)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MIKQRTLKRIVQATGVGLHTGKKVTLTLRPAPANTGVIYRRTDLNPPVDFPADA KSVRDTMLCTCLVNEHDVRISTVEHLNAALAGLGIDNIIVEVDAPEIPIMDGSAA PFVYLLLDAGIDELNCAKKFVRIKETVRVEDGDKWAEFKPYNGFSLDFTIDFNH PAIDASTQRYTLNFSADAFMRQISRARTFGFMRDIEYLQSRGLCLGGSFDCAIV VDDYRVLNEDGLRFEDEFVRHKMLDAIGDLFMCGHNIIGAFTAYKSGHALNNK LLQAVLAKQEAWYVTFEDDAKLPMAFRAPSMLVA
Research Area	Others
Source	E.coli
Target Names	IpxC
Expression Region	1-305aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	38.1 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.