




Recombinant Mouse Ras-related protein Rab-10 (Rab10)

Product Code	CSB-EP019152MOa0
Abbreviation	Recombinant Mouse Rab10 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61027
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAKKTYDLLFKLLLLIGDSGVGKTCVLFRRFSDDAFNTTFFISTIGIDFKIKTVELQGK KIKLQIWDTAGQERFHTITTSYYRGAMGIMLVYDITNGKSFENISKWLRNIDEHA NEDVERMLLGNKCDMDDKRVVPGKGEQIAREHGIRFFETSAKANINIEKAFLT LAEDILRKTPVKEPNSENVDISSGGGVTGWKSKCC
Research Area	Signal Transduction
Source	E.coli
Target Names	Rab10
Expression Region	1-200aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	26.6 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.