




Recombinant Human Atrial natriuretic peptide receptor 3 (NPR3), partial

Product Code	CSB-YP016025HU1f3
Abbreviation	Recombinant Human NPR3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17342
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GGVGGGGGGAGIGGGRQERREALPPQKIEVLVLLPQDDSYLFSLTRVRPAIEYA LRSVEGNGTGRLLPPGTRFQVAYEDSDCGNRALFSLVDRVAAARGAKPDLIL GPVCEYAAAPVARLASHWDLPLMSAGALAAGFQHKDSEYSHLTRVAPAYAKM GEMMLALFRHHHWSRAALVYSDDKLERNICYFTLEGVHEVFQEEGLHTSIYSF DETKDLLEDIVRNIQASERVVIMCASSDTIRSIMLVHRHGMTSGDYAFFNIEL FNSSSYGDGSWKRGDKHDFEAKQAYSSLQTVTLRTVKPEFEKFSMEVKSSV EKQGLNMEDYVNMFVEGFHDAILLYVLALHEVLRAGYSKKDGGKIIQQTWNRT FEGIAGQVSIANGDRYGD FSVIAMTDVEAGTQEVI G DYFGKEGRFEMRPNVK YPWGPLKLRIDENRIVEHTNSSPCKSSGGLEE
Research Area	Signal Transduction
Source	Yeast
Target Names	NPR3
Expression Region	27-481aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-Myc-tagged
Mol. Weight	53.4
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.