




Recombinant Bovine viral diarrhea virus Genome polyprotein, partial

Product Code	CSB-YP312464BKX1
Abbreviation	Recombinant Bovine viral diarrhea virus Genome polyprotein protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q01499
Product Type	Recombinant Proteins
Immunogen Species	Bovine viral diarrhea virus (strain SD-1) (BVDV) (Mucosal disease virus)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MELITNELLYKTYKQKPVGVEEPVYDQAGNPLFGERGAIHPQSTLKLPHKRGE RNVPTSLASLPKRGDCRSGNSKGPVSGIYLKPGPLFYQDYKGPVYHRAPLELF EEGSMCETTKRIGRVTGSDGKLYHIYICIDGCITVKSATRSHQRVLRWVHNRLD CPLWVTSC
Research Area	Others
Source	Yeast
Target Names	N/A
Protein Names	Recommended name: Genome polyprotein Cleaved into the following 13 chains: 1. N-terminal protease Short name= 2. N-pro EC= 3. 3.4.22.- Alternative name(s): Autoprotease p20 Capsid protein C E(rns) glycoprotein Alternative
Expression Region	1-168aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.4
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.