




Recombinant Hesperocyparis arizonica Pectate lyase 1

Product Code	CSB-YP871029COAD
Abbreviation	Recombinant Hesperocyparis arizonica Pectate lyase 1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCG9
Product Type	Recombinant Proteins
Immunogen Species	Hesperocyparis arizonica (Arizona cypress) (Cupressus arizonica)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DNPIDSCWRGDSNWDQNRMKLADCVVGFGSSTMGGKGGEIYTVTSSSEDPV NPTPGTLRYGATREKALWIIFSQNMNIKLQMPLYVAGYKTIDGRGADVHLGNG GPCLFMRKASHVILHGLHIHGCNTSVLGDVLVSESIGVEPVHAQDGDAITMRN VTNAWIDHNSLSDCSDGLIDVTLGSTGITISNNHFFNHHKVMMLLGHDDTYDDDK SMKVTVAFNQFGPNAGQRMPRARYGLVHVANNNYDQWNIYAIGGSSNPTILS EGNSFTAPNESYKKEVTKRIGCETTSACANWVWRSTRDAFTNGAYFVSSGKA EDTNIYNSNEAFKVENGNAAPQLTQNAGVVA
Research Area	Others
Source	Yeast
Target Names	N/A
Protein Names	Recommended name: Major pollen allergen Cup a 1 Alternative name(s): Allergen= Cup a 1
Expression Region	22-367aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	40.1
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.