




# Recombinant Vaccinia virus Entry-fusion complex associated protein OPG095 (OPG099), partial

<b>Product Code</b>	CSB-YP362178VAI1
<b>Abbreviation</b>	Recombinant Vaccinia virus VACWR088 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P07612
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	GAAASIQTTVNTLSERISSKLEQEANASQAQTKCDIEIGNFYIRQNHGCNLTVKN MCSADADAQLDAVLSAATETYSGLTPEQKAYVPAMFTAALNIQTSVNTVVRDF ENYVKQTCNSSAVVDNKLKIQNVIIDECYGAPGSPTNLEFINTGSSKGNCAIKAL MQLTTKATTQIAPKQVAGTG
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	VACWR088
<b>Protein Names</b>	Recommended name: Protein L1 Alternative name(s): Virion membrane protein M25
<b>Expression Region</b>	2-183aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	21.3 kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.