



Recombinant Human Putative teratocarcinoma-derived growth factor 3 (TDGF1P3)

Product Code	CSB-BP344744HUc7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51864
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	LGHQEFARPSRGDLAFRDDSIWPQEPAIRPRSSQRVLPNGIQHSKELNRTC CLNGGTCMLESFACPPSFYGRNCEHDVRKENCOSVPHDTWLPKKCSLCKC WHGQLRCFPQAFLPGCDGLVMDEHLVASRTPPELPPSARTTTFMLAGICLSIQS YY
Research Area	Others
Source	Baculovirus
Target Names	TDGF1P3
Expression Region	31-188aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	18.8
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.