




Recombinant Rotavirus A Intermediate capsid protein VP6

Product Code	CSB-BP365972RIX
Abbreviation	Recombinant Rotavirus A Intermediate capsid protein VP6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03530
Product Type	Recombinant Proteins
Immunogen Species	Rotavirus A (strain RVA/Human/United States/Wa/1974/G1P1A[8]) (RV-A)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEVLYSLSKTLKDARDKIVEGTLYSNVSDLIQQFNQMIVTMNGNDFQTGGIGNL PVRNWTDFDGLLGTTLLNLDANYVENARTIIEYFIDFIDNVCMDEMARESQRNG VAPQSEALRKLKAGIKFKRINFDNSSEYIENWNLQNRRTGFVFHKNIFPYSA SFTLNRSQPMHDNLMGTMWLNAGSEIQVAGFDYSCAINAPANIQQFEHIVQLR RALTTATITLLPDAERFSFPRVINSADGATTWFFNPVILRPNNVEVEFLLNGQIIN TYQARFGTIIARNFDAILLLFQLMRPPNMTPAVNALFPQAQPFQHHATVGLTLR IESAVCESVLADANETLLANVTAVRQEYAIPVGPVFPPGMNWTEELITNYSRE DNLQRVFTVASIRSMLIK
Research Area	Others
Source	Baculovirus
Target Names	N/A
Protein Names	Recommended name: Intermediate capsid protein VP6
Expression Region	1-397aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	48.8
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.