




Recombinant Anopheles gambiae Acetylcholinesterase (Ace), partial

Product Code	CSB-BP2916c7
Abbreviation	Recombinant Anopheles gambiae Acetylcholinesterase (Ace) protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q869C3
Product Type	Recombinant Proteins
Immunogen Species	Anopheles gambiae (African malaria mosquito)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	DNDPLVVNTDKGRIRGITVDAPSGKKVDVWLGIPIYAQPPVGPLRFRHPRPAEK WTGVLNTTTPPNSCVQIVDTVFGDFPGATMWNPNTPLSEDCLYINVVAPRPRP KNAAVMLWIFGGGFYSGTATLDVYDHRALASEENVIVVSLQYRVASLGFLFLG TPEAPGNAGLFDQNLALRWVRDNIHRFGGDP SRVT LFGESAGAVSVSLHLLS ALSRDLFQRAILQSGSPTAPWALVSREEATLRALRLAEAVGCPHEPSKLSDAV ECLRGKDPHVLVNNEWGTLGICEFPFVPVVDGAFLDETPQRSLASGRFKKTEI LTGSNTEEGYYFIIYYL TELLRKEEGVTVTREEFLQAVRELNPYVNGAARQAIVF EYTDWTEPDNPNSNRDALDKMVG DYHFTCNVNEFAQRYAEEGN NVMYLYT HRSKGNPWPRWTGVMHGDEIN YVFGEPLNPTLGYTEDEKDFSRKIMRYWSN FAKTGNPNPNTASSEFPEWPKHTAHGRHYLELGLNTSFVGRGPRLRQCAFW KKYLPQLVAATSNLPGPAPPSEPCCESS
Research Area	Others
Source	Baculovirus
Target Names	Ace
Expression Region	162-713aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	64.2
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.