




# Recombinant Human Kelch-like protein 12 (KLHL12)

<b>Product Code</b>	CSB-BP700650HU
<b>Abbreviation</b>	Recombinant Human KLHL12 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q53G59
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MGGIMAPKDIMTNTHAKSILNSMNSLRKSNTLCDVTLRVEQKDFPAHRIVLAAC SDYFCAMFTSELSEKGPYVDIQGLTASTMEILLDFVYTETVHVTVENVQELLP AACLLQLKGVKQACCEFLESQLDPSNCLGIRDFAETHNCVDLMQAAEVFSQK HFPEVVQHEEFILLSQGEVEKLIKDEIQVDSEEPVFEAVINWVKHAKKEREEES LPNLLQYVRMPLLTPRYITDVIDAEPFIRCSLQCRDLVDEAKKFHLRPELRSQM QGPRTRARLGANEVLLVGGFGSQSPIDVVEKYDPKTQEWSFLPSITRKRR YVASVSLHDRIYVIGGYDGRSRLSSVECLDYTADEGJVWYSVAPMNVRRGLA GATTLGDMIYVSGGFDGSRRTSMERYDPNIDQWSMLGDMQTAREGAGLVV ASGVIYCLGGYDGLNILNSVEKYDPHTGHWNTVTPMATKRSGAGVALLNDHIY VGGFDGTAHLSSVEAYNIRTDSWTTVTSMTPRCYVGATVLRGRLYAIAGYD GNSLLSSIECYDPIIDSWEVVTSMGTQRCDAGVCVLREK
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	Baculovirus
<b>Target Names</b>	KLHL12
<b>Protein Names</b>	Recommended name: Kelch-like protein 12 Alternative name(s): CUL3-interacting protein 1 DKIR homolog Short name= hDKIR
<b>Expression Region</b>	1-568aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	67.3
<b>Protein Length</b>	Full length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.