




Recombinant Mouse Neuronal pentraxin receptor (Nptxr), partial

Product Code	CSB-BP016031MO
Abbreviation	Recombinant Mouse Nptxr protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99J85
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ASVPLAASPARALPGGTDNASAASAAGGSGPQRSLHSAGGSAGPSVLPG EPAASVFPPPPVPLLSRFLCTPLAAACPSGAEQGDAAGERAELLLLQSTAEQL RQTALQQEARIRADRDRTIRELTGKLGRCESGLPRGLQDAGPRRDTMADGAWD SPALLLELEDAVRALRDRIERIEQELPARGNLSSAPAPAMPTALHSKMDELECQ LLAKVLALEKERAALSHGSHQQRQEVEKELNALQGRVAELEHGSSAYSPPDA FKVSIPIRNNYMYARVRKALPELYAFTACMCVRSRSGGSGQGTPFSYSVPGQ ANEIVLLEAGLEPMELLINDKVAQLPLSLKDSNWHHCISWTTTRDGLWSAYQDG ELRGSGENLAAWHPIKPHGILILGQEQDTLGGRFDATQAFVGDIAQFNLWDHA LTPAQVLGMANCTGPLMGNVLPWEDKLVEAFGGAKKAAFDVCKGRAKA
Research Area	Others
Source	Baculovirus
Target Names	Nptxr
Protein Names	Recommended name: Neuronal pentraxin receptor
Expression Region	24-493aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.8
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.