




Recombinant Mouse Extracellular tyrosine-protein kinase PKDCC (Pkdcc)

Product Code	CSB-BP735887MO
Abbreviation	Recombinant Mouse Pkdcc protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RJI4
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	PRPGQSPGSSAAPGPGRRGGRGELARQIRERYEEVQRYSRGGPGPGAGRP ERRRLMDLAPGGPGLQRPRPPRVRSPPDGAPGWPPAPGPGSPGPGPRLGC AALRNVSQAQYVGSYTKAVYRVRLPGGAVALKAVDFSGHDLGSCVREFGA RRGCYRLAAHKLLKEMVLLERLRHPNVLQLYGYCYQDSEGIPDTLTTITELGAP VEMIQLLQTSWEDRFRICLSLGRLLHHLAHSPLGSVTLDFRPRQFVLVNGELK VTDLDDARVEETPCTSSADCTLEFPARNFSLPCSAQGWCEGMNEKRNLYNAY RFFFTYLLPHSAPPSLRPLLDIVNATGELAWGVDETLAQLETALHLFRSGQYL QNSTSSRAEYQRIPDSAITQEDYRCWPSYHHGGCLLSVFNLAEAIDVCESHAQ CRAFVVTNQTWTGRKLVFFKTGWNQVVPDAGKTTYVKAPG
Research Area	Cell Biology
Source	Baculovirus
Target Names	Pkdcc
Protein Names	Recommended name: Protein kinase domain-containing protein, cytoplasmic EC= 2.7.-.- Alternative name(s): Protein kinase-like protein SgK493 Sugen kinase 493 Vertebrate lonesome kinase
Expression Region	33-492aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	54.4
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.