




Recombinant Mouse Enoyl-CoA hydratase, mitochondrial (Echs1)

Product Code	CSB-BP804339MO
Abbreviation	Recombinant Mouse Echs1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BH95
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ASGANFQYIITEKKGKNSSVGLIQLNRPKALNALCNGLIEELNQALETFEQDPAV GAIVLTGGDKAFAAGADIKEMQNRFTQDCYSSKFLSHWDHITRVKKPVIAAVN GYALGGGCELAMMCDIYAGEKAQFGQPEILLGTIPGAGGTQRLTRAVGKSLA MEMVLTGDRISAQDAKQAGLVSKIFPVEKLVVEEAIQCAEKIASNSKIVVAMAKE SVNAAFEMTLTEGNKLEKRLFYSTFATDDRREGMTAFVEKRKANFKDH
Research Area	Cancer
Source	Baculovirus
Target Names	Echs1
Protein Names	Recommended name: Enoyl-CoA hydratase, mitochondrial EC= 4.2.1.17 Alternative name(s): Enoyl-CoA hydratase 1 Short-chain enoyl-CoA hydratase Short name= SCEH
Expression Region	28-290aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	32.4
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.