




Recombinant Human Cytochrome P450 2C8 (CYP2C8)

Product Code	CSB-EP006418HU
Abbreviation	Recombinant Human CYP2C8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10632
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	RKLPPGPTPLPIIGNMLQIDVKDICKSFTNFSKVYGPVFTVYFGMNPVVFHGYE AVKEALIDNGEEFSGRGNPISQRITKGLGISSNGKRWKEIRRFSLTTLRNF GMRKSIDRVQEEAHCLVEELRKTASPCDPTFILGCAPCNVICSVVFQKRFDYK DQNFLTLMKRFNENFRILNSPWIVQVCNFPDLLIDCFPGTHNKVLKNVALTRSYI REKVKEHQASLDVNNPRDFIDCFLIKMEQEKDNQKSEFNENLVGTVADLFVA GTETTSTTLRYGLLLLLKHPEVTAKVQEEIDHVIGRHRSPCMQDRSHMPYTD VVHEIQRYSDLVPTGVPHAVTTDTKFRNYLIPKGTTIMALLTSVLHDDKEFPNP NIFDPGHFLDKNGNFKKSDYFMPFSAGKRICAGEGLARMELFLFLTTILQNFNL KSVDDLKLNLTAVTKGIVSLPPSYQICFIPV
Research Area	Cardiovascular
Source	E.coli
Target Names	CYP2C8
Protein Names	Recommended name: Cytochrome P450 2C8 EC= 1.14.14.1 Alternative name(s): CYP2C8 Cytochrome P450 IIC2 Cytochrome P450 MP-12 Cytochrome P450 MP-20 Cytochrome P450 form 1 S-mephenytoin 4-hydroxylase
Expression Region	27-490aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	60.1 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.