




Recombinant Mouse Transferrin receptor protein 1 (Tfrc), partial

Product Code	CSB-EP726792MO
Abbreviation	Recombinant Mouse Tfrc protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62351
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	CKRVEQKEECVKLAETEETDKSETMETEDVPTSSRLYWADLKLLSEKLNLSIEF ADTIKQLSQNTYTPREAGSQKDESLAYYIENQFHEFKFSKVWRDEHYVKIQVK SSIGQNMVTIVQSNGNLDPVESPEGYVAFSKPTEVSGKLVHANFGTKKDFEEL SYSVNGSLVIVRAGEITFAEKVANAQSFNAIGVLIYMDKNKFPVVEADLALFGHA HLGTGDPYTPGFPSFNHTQFPPSQSSGLPNIPVQTISRAAAELFGKMEGSCP ARWNIDSSCKLELSQNNVKNLIVKNVLRERRILNIFGVIKGYEEDRYVVVGAQ RDALGAGVAAKSSVGTGLLLKLAQVFSDMISKDGFRRSRSIIFASWTAGDFGA VGATEWLEGYLSSLHLKAFTYINLDKVVLTGTSNFKVSASPLLYTLMGKIMQDVK HPVDGKSLYRDSNWISKVEKLSFDNAAYPFLAYSGIPAVSFCFCEDADYPYLG TRLDTYEALTQKVPQLNQMVRTAAEVAGQLIKLTHDVELNLDYEMYSKLLSF MKDLNQFKTDIRDMGLSLQWLYSARGDYFRATSRLTDFHNAEKTNRVMREI NDRIMKVEYHFLSPYVSPRESPFRHIFWGSHTLSALVENLKLQRKNITAFNE TLFRNQLALATWTIQGVANALSGDIWNIDNEF
Research Area	Cancer
Source	E.coli
Target Names	Tfrc
Protein Names	Recommended name: Transferrin receptor protein 1 Short name= TR Short name= TfR Short name= TfR1 Short name= Trfr Alternative name(s): CD_antigen= CD71
Expression Region	89-763aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	80.0 kDa
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.