




Recombinant Rat E3 ubiquitin-protein ligase parkin (Park2)

Product Code	CSB-EP864348RAb1
Abbreviation	Recombinant Rat Park2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JK66
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MIVFVRFNSSYGFVPEVDS DTSIFQLKEVVAKRQGV PADQLRVIFAGKELQNH LTVQNC DLEQQSIVHIVQR PQRKSHETNASGGDKPQSTPEGSIWEPRSLTRVDL SSHILPADSVGLAVILD TDSKSDSEAARGPEAKPTYH SFFVYCKGPGCHKVQPG KLRVQCGTCRQATL TLAQGPSCWDDVLIPNRMSGECQSPDCPGTRA EFFFFKC GAHPTSDK DTSVALNLITNNSRSIPCI ACTDVRNPVLFQCNHRHVICLDCFHL YCVTRLNDRQFVHDAQLGYS LPCVAGCPNSLIKELH HFRILGEEQYNRYQQYG AEECVLQMGGVLCPRPGCGAGLLPEQGQKKVTCEGGNGLGCGFVFCRDCK EAYHEGECDSMFEASGATSQAYRVDQRAAEQARWEEASKETIKKTTKPCPRC NVPIEKNGGCMHMKCPQPQCKLEWCWNCGCEWNRACMGDHWFDV
Research Area	Cell Biology
Source	E.coli
Target Names	Prkn
Expression Region	1-465aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	56.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.