




Recombinant Rabies virus Glycoprotein (G), partial

Product Code	CSB-EP325676RAFe1
Abbreviation	Recombinant Rabies virus Glycoprotein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19462
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rabies virus (strain HEP-Flury) (RABV)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KFPIYTIPDKLGPWSPIDLHHLSCPNNLVVEDEGCTNLSGFSYMEKLVGYISAIAIK VNGFTCTGVVTEAETYTNFVGIVTTTTFKRKHFRPTPDACRAAYNWKMAGDPR YEESLHNPYPDYHWRRTVKTTKESLVIISPSVTDLDPYDKSLHSRVFPGGNCS GITVSSTYCSTNHDYTIWMPENLRLGTSCDIFTHSRGKRASKGDKTCGFVDER GLYKSLKGACKLKLKCGVLGLRLMDGTWVAMQTSDETKWCPPGQLVNLHDFR SDEIEHLVEEELVKKREECLDALESIMTTKSVSFRRLSHLRKLVPFGFKAYTIFN KTLMEADAHYKSVQWNEIIPSKGCLRVGERCHPHVNGVFFNGIILGSDGHVLI PEMQSSLLQQHMELLESSVIPLMHPLADPSTVFKDGDVEVDFVEVHLPDVHKQ VSGVDLGLPKWGKY
Research Area	Signal Transduction
Source	E.coli
Target Names	G
Expression Region	20-459aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	49.7 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.