




# Recombinant Mouse E3 ubiquitin-protein ligase parkin (Prkn)

<b>Product Code</b>	CSB-EP893911MO
<b>Abbreviation</b>	Recombinant Mouse Prkn protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9WVS6
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus(Mouse)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MIVFVRFNSSYGFVPEVDS DTSILQLKEVVAKRQGVPADQLRVIFAGKELPNHL TVQNC DLEQQSIVHIVQRPRRRSHETNASGGDEPQSTSEGSIWESRSLTRVDL SSHTLPVDSVGLAVILD TDSKRDSEAARGPVKPTYNSFFIYCKGPCHKVQPGK LRVQCGTCKQATLTLAQQPSCWDDVLIPNRMSGECQSPDCPGTRAEFFFKC GAHPTSDKDTSVALNLITSNRRSIPCIACTDVRSPVLVFQCNHRHVICLDCFHLY CVTRLNDRQFVHDAQLGYSLPCVAGCPNSLIKELHHFRILGEEQYTRYQQYGA EECVLQMGGVLCPRPGCGAGLLPEQQQRKVTCEGGNGLGCGFVFCRDCKE AYHEGDCDSLLEPSGATSQAYRVDKRAAEQARWEEASKETIKKTTKPCPRCN VPIEKNGGCMHMKCPQPQCKLEWCWNCGCEWNRACMGDHWFDV
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	Prkn
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase parkin EC= 6.3.2.-
<b>Expression Region</b>	1-464aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	79.1 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.