




Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2

Product Code	CSB-EP310896DNLa0
Abbreviation	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P92918
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Apium graveolens (Celery)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MGVQKTVEAPSTVSAEKMYQGFLDMDTVFPKVLPLQLIKSVEILEGDGGVGT VKLVHLGEATEYTTMKQKVDVIDKAGLAYTYTTIGGDILVDVLESVVNEFVVVPT DGGCIVKNTTIYNTKGDAVLPEDKIKEATEKSALAFKAVEAYLLANLQFLA
Research Area	Others
Source	E.coli
Target Names	N/A
Expression Region	1-159aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	23.1 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.