




Recombinant Human Solute carrier organic anion transporter family member 2B1 (SLCO2B1), partial

Product Code	CSB-YP021746HU
Abbreviation	Recombinant Human SLCO2B1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94956
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	FFIGCSSHQIAGITHQTSAHPGLELSPSCMEACSCPLDGFNPVCDPSTRVEYIT PCHAGCSSWVVQDALDNSQVFYTNCSCVVEGNPVLGAGSCDSTCSHLVVPF
Research Area	Others
Source	Yeast
Target Names	SLCO2B1
Protein Names	Recommended name: Solute carrier organic anion transporter family member 2B1 Alternative name(s): Organic anion transporter B Short name= OATP-B Organic anion transporter polypeptide-related protein 2 Short name= OATP-RP2 Sho
Expression Region	461-564aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	13
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.