




Recombinant Human Putative teratocarcinoma-derived growth factor 3 (TDGF1P3), partial

Product Code	CSB-YP344744HU
Abbreviation	Recombinant Human TDGF1P3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51864
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	LGHQEFARPSRGDLAFRDDSIWPQEEP AIRPRSSQRVLP MGIQHSKELNRTC CLNGGTCMLESFACPPSFYGRNCEHDVRKENCGSVPHDTWLPKCKSLCKC WHGQLRCFPQAFLPGCDGLVMDEHLVASRTPPELPPSARTTTFMLAGICLSIQS YY
Research Area	Others
Source	Yeast
Target Names	TDGF1P3
Protein Names	Recommended name: Putative teratocarcinoma-derived growth factor 3 Alternative name(s): Cripto-3 growth factor Epidermal growth factor-like cripto protein CR3 Teratocarcinoma-derived growth factor 1 pseudogene 3
Expression Region	31-188aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	18.6
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.