




Recombinant Human YjeF N-terminal domain-containing protein 3 (YJEFN3)

Product Code	CSB-YP026264HU
Abbreviation	Recombinant Human YJEFN3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	May play a role in spermiogenesis and oogenesis.
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSSAAGPDPSEAPEERHFLRALELQPPLADMGRAELSSNATTSLVQRRKQAW GRQSWLEQIWNAGPVCQSTAEAAALERELLEDYRFGRQQVLVELCGHASAVAV TKAFPLPALSRKQRTVLVVCGPENGAVGLVCAHRLRVFEYEPTIFYPTRSLDL LHRDLTTQCEKMDIPFLSYLPTEVQLINEAYGLVVDVAVLGPVGEVGGPCT RALATLKLLSIPLVSLDIPSGWDAETGSDSEDGLRPDVLVSLAAPKRCAGRFSG RHHFVAGRFPDDVRRKFALRLPGYTGTCVAAL
Research Area	Others
Source	Yeast
Target Names	YJEFN3
Protein Names	Recommended name: YjeF N-terminal domain-containing protein 3 Short name= YjeF_N3 Short name= hYjeF_N3
Expression Region	1-299aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	33.2
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.