




Recombinant Human Protein furry homolog-like (FRYL), partial

Product Code	CSB-YP009007HU
Abbreviation	Recombinant Human FRYL protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94915
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	KLLSRSLQRGEDLQFDQLISSMSSVAEHCLPSLLRTLFDWYRRQNGTEDES YEYRPRSSTKSKGDEQQRERDYLLERRDLAVDFIFCLVLVEVLKQIPVHPVPDPLV HEVLNLAFFKHFKHKEGYSGTNTGNVHIIADLYAEVIGVLAQSKFQAVRKKFVTE LKELRQKEQSPHVVSQSVISLIMGKFFRVKMYPVDFEASFQFMQECAQYFLE VKDKDIKHALAGLFVEILIPVAAAVKNEVNPCLKNFVEMLYQTTFELSSRKKHS LALYPLITCLLCVSQKQFFLNNWHIFLQNCLSHLKKNKDPKMSRVALESYRLLW VYVIRIKCESNTVTQSRMSIVSALFPKGSRSVVRDTPNLNIFVKIIQFIAQERLD FAMKEIIFDLLSVGKSTKTFTINPERMNIPLRVFLVIADSLQKQDGEPPMPTTGVI LPSGNTLRVKKIFLNKTLTDEEAKVIGMSVYYPQVRKALDSILRHLDKEVGRPM CMTSVQMSNKEPEDMITGERKPKIDLFRTCIAAIPRLPDGMSRTDLIELLARLTI HMDEELRALAFNTLQALMLDFPDWREDVLSGFVYFIVREVTDVHPTLLDNAVK MLVQLINQWKQAAQMHNKNQDTQHGVANGASHPPPLERSPYSNVHVVEGF ALVILCSSRPATRRLAVSVLREIRALFALLEIPKGDDELAIDVMDRLSP
Research Area	Others
Source	Yeast
Target Names	FRYL
Protein Names	Recommended name: Protein furry homolog-like Alternative name(s): ALL1-fused gene from chromosome 4p12 protein
Expression Region	42-739aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-Myc-tagged
Mol. Weight	83.40000000000001
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.