




# Recombinant Mouse Leucine-rich repeat LGI family member 3 (Lgi3)

<b>Product Code</b>	CSB-EP812997MO
<b>Abbreviation</b>	Recombinant Mouse Lgi3 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K406
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus(Mouse)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	KRPPKTPPCPPSCSCTRDTAFVCVDSKSVPKNLPSEVISLTLVNAAFSEIQDGAF SHLPLLQFLLLNSNKFTLIGDNAFIGLSHLQYLFIEENNDIWALSKFTFRGLKSLTH LSLANNNLQTLPRDIFRPLDILSDLDLRGNALNCDCVKVWLVEWLAHTNNTTVAP IYCASPPRFQEHKVQDLPLREFDCITTDVLYQTLSPFAVSAEPFLYSSDLYLAL AQPASACTILKWDYVERQLRDYDRIPAPSAVHCKPMVVDGQLYVVVAQLFG GSYIYHWDPNTRFTKLQDIDPQRVRKPNLDLEAFRIDGDWFFAVADSSKAGAT SLYRWHQNGFYSHQALHAWHRDTDLEFVDGEGKPRILIVSSSSQAPVIYQWSR SQKQFVAQGEVTQVPDAQAVKHFRAGRDSYLCLSRIGDSKILRWEGTRFSE VQALPSRGLALQPFLVGGHRYLALGSDFSFTQIQWDEGRQKFVRFQELAV QAPRAFICYMPAGDAQLLLAPSFKGQTLVYRHVVVDLSA
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	Lgi3
<b>Protein Names</b>	Recommended name: Leucine-rich repeat LGI family member 3 Alternative name(s): Leubrin Leucine-rich glioma-inactivated protein 3
<b>Expression Region</b>	31-548aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	62.7 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.