




Recombinant Mouse Chemokine-like protein TAFA-1 (Tafa1)

Product Code	CSB-EP763174MO
Abbreviation	Recombinant Mouse Tafa1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TPG8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MLLCHGSLQHTFQQHHLHRPEGGTCEVIAAHRCCNKNRIEERSQTVKCSCLP GKVAGTTRNRPSCVDASIVIGKWWCEMEPCLEGEECKTLPDNSGWMCATGN KIKTTRIHPRT
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Fam19a1
Protein Names	Recommended name: Protein FAM19A1Alternative name(s): Chemokine-like protein TAFA-1
Expression Region	20-133aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	20.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.