




Recombinant Human CD81 antigen (CD81), partial

Product Code	CSB-MP004960HU
Abbreviation	Recombinant Human CD81 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P60033
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	FVNKDQIAKDVKQFYDQALQQAVVDDANNAAKAVVKTFFHETLDCCGSSTLTAL TTSVLKNNLCPSPGSNIISNLFKEDCHQKIDDLFSGK
Research Area	Immunology
Source	Mammalian cell
Target Names	CD81
Protein Names	Recommended name: CD81 antigen Alternative name(s): 26 kDa cell surface protein TAPA-1 Target of the antiproliferative antibody 1 Tetraspanin-28 Short name= Tspan-28 CD_antigen= CD81
Expression Region	113-201aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	38.7
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.