




# Recombinant Human Tumor necrosis factor ligand superfamily member 14 (TNFSF14), partial

<b>Product Code</b>	CSB-MP023991HU
<b>Abbreviation</b>	Recombinant Human TNFSF14 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O43557
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	DGPAGSWEQLIQERRSHEVNPA AHLTGANSSLTGSGGPLLWETQLGLAFLRG LSYHDGALVVTKAGYYYIYSKVQLGGVGCPLGLASTITHGLYKRTPRYP EEELEL LVSQQSPCGRATSSSRVWWDSSFLGGVVHLEAGEKVVVVRV LDERLVRLRDG TRSYFGAFMV
<b>Research Area</b>	Immunology
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TNFSF14
<b>Expression Region</b>	74-240aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal hFc1-Myc-tagged
<b>Mol. Weight</b>	48.3
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.