




Recombinant Human Keratin, type II cytoskeletal 8 (KRT8)

Product Code	CSB-BP012574HU(A4)
Abbreviation	Recombinant Human KRT8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05787
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSIRVTQKSYKVSTSGPRAFSSRSYTSGPGSRISSSSF SRVGSSNFRGGLGG GYGGASGMGGITAVTVNQSLLSPLVLEVDPNIAVVRTQEKEQIKTLNNKFASFI DKVRFLEQQNKMLETKWSLLQQKTARSNMDNMFESYINNLRRLQLETLGQEK LKLEAELGNMQGLVEDFKNKYEDEINKRTEMENEFVLIKKDVDEAYMNKVELE SRLEGLTDEINFLRQLYEEIRELQSQISDTSVVLSDNRSRLDMSIIAEVKAQ YEDIANRSRAEAESMYQIKYEELQSLAGKHGDDLRRRTKTEISEMNRNISRQLAE IEGLKGQRASLEAAIADAEQRGELAIKDANAKLSELEAALQRAKQDMARQLRE YQELMNVKLALDIEIATYRKLLERGEESRLESGMQNMSIHTKTTSGYAGGLSSAY GGLTSPGLSYSLSGSSFGSGAGSSSFRTSSSRVVVKIETRDGKLVSESSDV LPK
Research Area	Signal Transduction
Source	Baculovirus
Target Names	KRT8
Expression Region	1-483aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	57.6
Protein Length	Full length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.