




Recombinant Human Adiponectin receptor protein 1 (ADIPOR1), partial

Product Code	CSB-CF001367HU2
Abbreviation	Recombinant Human ADIPOR1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96A54
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EGRWRVIPYDVLDPDWLKDNDYLLHGHHRPPMPSFRACFKSIFRIHTETGNIWTH LLGFVLFLLFLGILTMLRPNMYFMAPLQEKVVFGMFFLGAVLCLSFSLFHTVYC HSEKVSRTFSKLDYSGIALLIMGSFVPWLYYSFYCSPQPRLIYLSIVCVLGISAIIV AQWDRFATPKHRQTRAGVFLGLGLSGVVPTMHFTIAEGFVKATTVGQMGWF FLMAVMIYTGAGLYAARIPERFFPGKFDIWFQSHQIFHVLVVAAAFVHFYGVSN LQEFRYGLEGGCTDDTLL
Research Area	Cardiovascular
Source	in vitro E.coli expression system
Target Names	ADIPOR1
Expression Region	89-375aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Flag-tagged
Mol. Weight	34.8 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.