



Recombinant Tick-borne encephalitis virus

Genome polyprotein, partial

Product Code	CSB-EP3056GKR
Abbreviation	Recombinant Tick-borne encephalitis virus Genome polyprotein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A075TBC5
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Tick-borne encephalitis virus
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	GLTYTMCDKTKFTWKRIPTDSGHDTVMEVAFSGTKPCRIPVRAVAHGSPDV NVAMLITP NPTIETNGGGFIEMQLPPGDNIIYVGELSHQWFQKGSSI
Research Area	Others
Source	E.coli
Expression Region	302-400aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	26.1 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.