




Recombinant Human AT-rich interactive domain-containing protein 4B (ARID4B), partial

Product Code	CSB-EP702965HU
Abbreviation	Recombinant Human ARID4B protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4LE39
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KQIDELLGKVVCVDYISLDKKKALWFPALVPCDCSDEIAVKKDNILVRSFKDG KFTSVPRKDVHEITSDTAPKPDVAVLKQAFEQALEFHKSRTIPANWKTELKEDSS SSEAEEDDEKEKEDNSSEEEIEPFPEERENFLQQLYKFMEDRGTP NKRPVLGYRNLNLFKLFRLVHKLGGFDNIESGAVWKQVYQDLGIPVLNSAAGY NVKAYKKYLYGFEEYCRSANIEFQMALPEKVVNK
Research Area	Others
Source	E.coli
Target Names	ARID4B
Protein Names	Recommended name: AT-rich interactive domain-containing protein 4B Short name= ARID domain-containing protein 4B Alternative name(s): 180 kDa Sin3-associated polypeptide Short name= Sin3-associated polypeptide p180 Breast cancer-associ
Expression Region	167-415aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	32.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.