




# Recombinant Mouse Zinc finger protein Helios (Ikzf2)

<b>Product Code</b>	CSB-EP011576MO
<b>Abbreviation</b>	Recombinant Mouse Ikzf2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P81183
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus(Mouse)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	METDAIDGYITCDNELSPEGEHANMAIDLTSSTPNGQHASPSHMTSTNSVKLE MQSDEECDRQPLSREDEIRGHDEGSSLEEPLIESSEVADNRKVQDLQGEggi RLPNGKLCQDVCGMVCIGPNVLMVHKRSHTGERPFHCNQCASFTQKGNLL RHIKLSGEEKPFKCPFCYACRRRDALTGHLRTHSVGKPHKCNYCGRSYKQR SSLEEHKERCHNYLQNVSMEEAAGQVMSSHVPPMEDCKEQEPIMDNNISLVPF ERPAVIEKLTANMGKRKSSTPQKFVGEKLMRFSYPDIHFDMNLTYEKEAELMQ SHMMDQAINNAITYLGAEALHPLMQHAPSTIAEVAPVISSAYSQVYHPNRIERPI SRETSDSHENMDGPISLIRPKSRPQEREASPSNSCLDSTDSESSHDDRQSY QGNPALNPKRKQSPAYMKEDVKALDATKAPKGS�KDIYKVFNGEGEQIRAFKC EHCrvLFLDHVMYTIHMGCHGYRDPLECNICGYRSQDRYEFSSHIVRGEHTFH
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	E.coli
<b>Target Names</b>	Ikzf2
<b>Protein Names</b>	Recommended name: Zinc finger protein Helios Alternative name(s): Ikaros family zinc finger protein 2
<b>Expression Region</b>	1-526aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	66.8 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.