




# Recombinant Acinetobacter baumannii UDP-3-O-acyl-N-acetylglucosamine deacetylase (IpxC)

<b>Product Code</b>	CSB-EP460177AWOe1
<b>Abbreviation</b>	Recombinant Acinetobacter baumannii IpxC protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B2I1I9
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Acinetobacter baumannii(strainACICU)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MVKQRTLNRVVKASGIGLHSGQKVMINFIPTVDGGIVFRRIDLDPVDIPANAL LIQEAFCMCSNLVTGDIKVGTEIEHVMSAIAAGLGIDNLIVEVSASEVPIMDGSAGPFI YLLMQGGLREQDAPKKFIKILKPVEALIDDKKAIFSPHNGFQLNFTIDFDHPAFA KEYQSATIDFSTETFFVYEVSEARTFGFMKDLDYLNKANNLALGASLDNAIGVDDT GVVNEEGLRFADEFVRHKILDAVGDLYLLGHQIIAKFDGYKSGHALNNQLLRNV QSDPSNYEIVTFDDEKDCPIPYVSVT
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	IpxC
<b>Expression Region</b>	1-300aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	33.1 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.