




# Recombinant Human Keratin, type I cytoskeletal 18 (KRT18)

<b>Product Code</b>	CSB-EP012519HUe1
<b>Abbreviation</b>	Recombinant Human KRT18 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P05783
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	SFTTRSTFSTNYRSLGSVQAPSYGARPVSSAASVYAGAGGSGSRISVSRSTSF RGGMGSGGLATGIAGGLAGMGGIQNEKETMQSLNDRLLASYLDRVRSLETENR RLESKIREHLEKKGQVRDWSHYFKIIEDLRAQIFANTVDNARIVLQIDNARLAA DDFRVKYETELAMRQSVENDIHGLRKVIDDTNITRLQLETEIEALKEELLFMKKN HEEEVKGLQAQIASSGLTVEVDAPKSQDLAKIMADIRAQYDELARKNREELDK YWSQQIEESTTVTTQSAEVGAAETTLTELRRTVQSLEIDLDSMRNLKASLENS LREVEARYALQMEQLNGILLHLESELAQTRAEGQRQAQEQYEALLNIKVKLEAEI ATYRRLLEDGEDFNLGDALDSSNSMQTIQKTTTTRIVDGKVVSETNDTKVLRH
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	KRT18
<b>Expression Region</b>	2-430aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	48.1 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.