




Recombinant Mouse Beta-tectorin (Tectb)

Product Code	CSB-EP023371MOB9
Abbreviation	Recombinant Mouse Tectb protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08524
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KSCTPNKADVILVFCYPKTIITKIPECPYGWEVHQLALGGGLCYNGVHEGGYYQF VIPDLSPKNKSYCGTQSEYKPPIYHFYSHIVSNDSTVIVKNQPVNYSFSCTYHS TYLVNQAAFDQRVATVHVKNKSGMTFESQLSLNFYTNAKFSTKKEAPFVLETS EIGSDLFAGVEAKGLSVRFKVVNLNSCWATPSADFMYPQLQWQLINKGCPTDET LVHENGKDHRATFQFNAFRFQNIPLSKVWLHCETFICDSEKLSCPVNCDKRRK RMLRDQTGGVLVVELSLRSRA
Research Area	Others
Source	E.coli
Target Names	Tectb
Expression Region	18-305aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-B2M-JD-tagged and C-terminal Myc-tagged
Mol. Weight	39.4 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.