




# Recombinant Human Laminin subunit gamma-2 (LAMC2), partial

<b>Product Code</b>	CSB-EP618803HU
<b>Abbreviation</b>	Recombinant Human LAMC2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13753
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	NCQGGGACDPDTGDCYSGDENPDIECADCPIGFYNDPHDPRSCCKPCPCHNG FSCSVMPETEEVVCNNCPPGVTGARCELCADGYFGDPFGEHGPVRPCQPCQ CNNNVDPASGNCDRLTGRCLKCIHNTAGIYCDQCKAGYFGDPLAPNPADKC RACNCNPMGSEPVGCRSD
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	LAMC2
<b>Protein Names</b>	Recommended name: Laminin subunit gamma-2 Alternative name(s): Cell-scattering factor 140 kDa subunit Short name= CSF 140 kDa subunit Epiligrin subunit gamma Kalinin subunit gamma Kalinin/nicein/epiligrin 100 kDa subunit Lad
<b>Expression Region</b>	417-588aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-GST-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	53.3 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.