




Recombinant Triticum aestivum Glutenin, high molecular weight subunit DX5 (GLU-1D-1D), partial

Product Code	CSB-EP319188TQN
Abbreviation	Recombinant Triticum aestivum GLU-1D-1D protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10388
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Triticum aestivum(Wheat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EGEASEQLQCERELQELQERELKACQQVMDQQLRDISPECHPVVWSPVAGQ YEQQIVVPPKGGSFYPGETTPPQQLQQRIFWGIPALLKRYYPSTCPQQVSY PGQASPRPGQGQPGQGQGGYYPTSPQQPGQWQQPEQGQPRYYPTSP QQSGQLQQPAQGQPGQGQGG
Research Area	Others
Source	E.coli
Target Names	GLU-1D-1D
Protein Names	Recommended name: Glutenin, high molecular weight subunit DX5
Expression Region	22-194aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	20.8 kDa
Protein Length	partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.