




Recombinant Mycoplasma pneumoniae P30 adhesin (p30), partial

Product Code	CSB-EP301931MLW1
Abbreviation	Recombinant Mycoplasma pneumoniae P30 adhesin protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P75330
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	RLLEEKERQEQLAEQLQRISAQQEEQQALEQQAAAEAHAEAEVEPAPQPVPV PPQPQVQINFGPRTGFPPQPGMAPRPGMPHPGMAPRPGFPPQPGMAPRP GMPPHPGMAPRPGFPPQPGMAPRPGMPHPGMAPRPGFPPQPGMAPRPG MQPPRPGMPQPGFPPKR
Research Area	Neuroscience
Source	E.coli
Target Names	p30
Protein Names	Recommended name: P30 adhesin Alternative name(s): 30 kDa adhesin-related protein Cytadhesin P30
Expression Region	106-274aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-GST-tagged
Mol. Weight	49.6 kDa
Protein Length	partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.