



# Recombinant Human Runt-related transcription factor 1 (Runx1)

<b>Product Code</b>	CSB-EP020592HUb1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q01196
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MRIPVDASTSRRFTPPSTALSPGKMSEALPLGAPDAGAALAGKLRSGDRSMV EVLADHPGELVRTDSPNFLCSVLPTHWRCNKTLPIAFKVVVALGDVPDGLTVTV MAGNDENYSAELRNATAAMKNQVARFNDLRFVGRSGRGKSFTLTITVFTNPP QVATYHRAIKITVDGPREPRRHRQKLLDDQTKPGSLSFSERLSELEQLRRTAMR VSPHHPAPTPNPRASLNHSTAFNPQPQSQMQDTRQIQSPPPWSYDQSYQYL GSIASPSVHPATPISPGRASGMTTLSAELSSRLSTAPDLTAFSDPRQFPALPSIS DPRMHYPGAFTYSPTPVTSGIGIGMSAMGSATRYHTYLPPPYPGSSQAQGGP FQASSPSYHLYYGASAGSYQFSMVGGERSPPRILPPCTNASTGSALLNPPLPN QSDVVEAEGSHSNSPTNMAPSARLEEAVWRPY
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	Runx1
<b>Expression Region</b>	1-453aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	53.8 kDa
<b>Protein Length</b>	Full Length
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.