




Recombinant Human Filamin-A-interacting protein 1 (FILIP1), partial

Product Code	CSB-EP800226HU
Abbreviation	Recombinant Human FILIP1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z7B0
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KSLELSKEDLIQLLSIMEGELQAREDVIHMLKTEKTKPEVLEAHYGSAEPEKVLR VLHRDAILAQEKSIGEDVYEKPISELDRLEEKQKETYRRMLEQLLLAEKCHRR T VYELENEKHKHTDYMNKSDDFTNLLEQERERLKKLLEQEKEYQARKEKENAK RLNKLRLDELVKLKSFALMLVDERQMHIEQ
Research Area	Others
Source	E.coli
Target Names	FILIP1
Protein Names	Recommended name: Filamin-A-interacting protein 1 Short name= FILIP
Expression Region	71-260aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	30.2 kDa
Protein Length	partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.