




Recombinant Human Thrombomodulin (THBD), partial

Product Code	CSB-EP023486HU2
Abbreviation	Recombinant Human THBD protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07204
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	RGADFQALPVGSSAAVAPLGLQLMCTAPPGAVQGHWAREAPGAWDCSVEN GGCEHACNAIPGAPRCQCPAGAALQADGRSCTASATQSCNDLCEHFCVNP DQPGSYSCMCETGYRLAADQHRCEDVDDCILEPSPCPQRCVNTQGGFECHC YPNYDLVDGECVEPVDPCFRANCEYQCQPLNQTSYLCVCAEGFAPIPHEPHR CQMFCNQATACPADCDPNTQASCECPEGYILDDGFICTDIDECENGGFCSGVC HNLPGTFECICGPDSALARHIGTDCCDSGKVDGGDSGSGEPPPSPT
Research Area	Cardiovascular
Source	E.coli
Target Names	THBD
Expression Region	200-500aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	47.0 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.