




Recombinant Mouse Heterogeneous nuclear ribonucleoproteins A2/B1 (Hnrnpa2b1)

Product Code	CSB-EP010602MOb0
Abbreviation	Recombinant Mouse Hnrnpa2b1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O88569
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEKTLETVPLERKKREKEQFRKLFIGGLSFETTEESLRNYYEQWGKLTDCVVM RDPASKRSRGFGFVTFSSMAEVDAAAMAARPHSIDGRVVEPKRAVAREESGKP GAHVTVKKLFVGGIKEDTEEHHLRDYFEEYGKIDTIEITDRQSGKKRGFGFVTF DDHDPVDKIVLQKYHTINGHNAEVRKALSRQEMQEVQSSRSRGGNFGFGD SRGGGNFGPGPGSNFRGGSDGYGSGRGFGDGYNGYGGGPGGGNFGGS PGYGGGRGGYGGGGPGYGNQGGGYGGGYDNYGGGNYGSGSYNDFGNYN QQPSNYGPMKSGNFGGSRNMGGPYGGGNYGPGGSGGSGGYGGRSRY
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Hnrnpa2b1
Expression Region	1-353aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	43.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.